

# Connectivity

## Linking to other devices and technologies

Connect your hearing devices to other technologies and accessories.

### Examples of what this can help you with:

- Stream music to my hearing aids from my phone
- Be able to hear the TV
- Talk on the telephone effectively

## Telecoil



A small coil inside the hearing aid that picks up sounds directly from compatible phones and hearing loops. These may be used in public places, such as concert halls and theatres.

## Mobile Phone and Device Connectivity



Connect your hearing aids with Bluetooth® devices. May be used for phone calls, audio streaming and hearing aid control.

## Accessory Connectivity



Links your hearing aids to hearing aid accessories. This can include remote microphones and TV streamers.

## Phone Assistive Technology



Detects when a phone is being used and adjusts hearing aid settings for you.



Scan for more detailed technical information



v1.2 2/11/23